non sterile

Revision: 241005

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INSTRUCTIONS FOR USE 4. (GB) Warning: The implants described may only be used in combination with erbrich products!

IMPLANT – HIP ENDOPROSTHESIS non-sterile (HEAD PROSTHESES OF THE MOORE AND THOMPSON BRAND)

IMPLANT MATERIAL:
Implants made from steel in accordance with ISO 5832 – 1 :1997 Stainless steel manufactured by erbrich-Instrumente GmbH.

The surface of these implants is chemically passive: not magnetic. These implants can be combined with the standardised material provided that its composition lies within the analysis stipulated in the ISO 5832 – 1:1997 standard and the required specifications. Material alloys with chrome-nickel (implant steel) components can trigger a so-called chrome-nickel allergy by their nature; there are no other known biocompatible hazards!

We would point out that implants only function correctly if the following basic rules are observed:

-When selecting implants (prostheses),care must be taken to choose the appropriate implant (hip prosthesis) on the basis of the weight and level of activity of the patient, and the hip joint to be treated.
-Note that by making the correct choice of biomechanics the forces to be transferred by the implants remain low. Damaged or deformed implants must not be implanted.
-The re-use of implants is absolutely forbidden (single use).
-We strongly advise you to inform patients of the advantages and disadvantages of implants.
-Excessive strain as a result of body weight and the level of activity of the patient should be avoided due to the limited strength of the implants. Failure to observe these precautionary measures can have serious consequences for the healing process.

PREOPERATIVE PLANNING
Preoperative planning provides important information on the suitable prosthesis. Different sizes of prosthesis should be available in case other sizes are required or the planned prosthesis cannot be used.

INDICATIONS
Advanced attrition of the hip joint owing to degenerative, post-traumatic or rheumatic arthritis; fracture or avascular necrosis of the femoral head; resulting condition of earlier operations e.g. osteosynthesis, joint reconstruction, arthrodesis.

CONTRAINDICATIONS
Acute and chronic infections; muscle, nerve or vascular diseases that endanger the affected extremity; a lack of bony tissue or poor bone quality (e.g. severe osteoporosis); local bone tumours. Systemic diseases and metabolic disorders; infections and falls; drug dependence; obesity; highly physical activities together with activities involving extreme vibrations which can lead to overtraining of the implant (prosthesis).

NOTES
The user should record and keep all information provided to the patient. It should be checked before use whether the patient tolerates the material to be implanted. The implants described in these instructions for use may only be used (implanted) by surgeons with the appropriate experience.

POSSIBLE NEGATIVE EFFECTS (RISKS)
Displacement and loosening of the prosthesis; luxation of the prosthesis; infection; venous thrombosis and pulmonary embolism; cardiovascular disorders; haematomas

CLEANING AND STERILISATION OF UNSTERILE IMPLANTS
Before use, the implants must be sterilised using a validated procedure (autoclaved).

ERBRICH – INSTRUMENTE GMBH VALIDATED PROCEDURES
STERILISER: Steam autoclave: Temperature = 121° to 123° Celsius; Pressure = 1 – 1.2 bar (15 to 17 psi) with an exposure time of at least 30 minutes in the packaged state.
STERILISER: Autoclave with prevacuum: During the prevacuum phase the air is extracted from the chamber before the steam flows in. Standard circulation for packaged objects: 132° to 135° Celsius; Pressure = 2 – 3 bar (27 to 30 psi) with an exposure time of at least 4 to 10 minutes.

IMPORTANT NOTES: Every time the implants are used or sterilised you must first check that they function perfectly!
In the event of failure to observe these instructions or demonstrable violation, no guarantee or compensation for damages can be provided.

erbrich - instrumente GmbH
Eugenstr. 33 · D-78532 Tuttlingen
Tel.: +49 74 61 / 96 57 30
Fax: +49 74 61 / 9657 320
E-mail: info@erbrich-instruments.com
Website: www.erbrich-instruments.com

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Impianti in acciaio ai sensi di ISO 5832 – 1:1997 Acciaio inossidabile prodotto dalla ditta erbrich-Instrumente GmbH.